

## Product Datasheet

### Rabbit Anti-Human IgG F(ab')<sub>2</sub> - AP (orb1878877)

<b>Description</b>	Rabbit Polyclonal Secondary Antibody to Human IgG F(ab') <sub>2</sub> AP labeled
<b>Reactivity</b>	Human
<b>Conjugation</b>	AP
<b>Tested Applications</b>	ELISA, IH, WB
<b>Immunogen</b>	Human IgG F(ab') <sub>2</sub>
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Dilution Range</b>	E (1/10000 - 1/50000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
<b>Expiration Date</b>	12 months from date of receipt.