

Product Datasheet

Goat Anti-Rat IgG (H&L) - FITC (orb1878923)

Description	Goat Polyclonal Secondary Antibody to Rat IgG (H&L) FITC
Reactivity	Rat
Conjugation	FITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Rat IgG
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	IF (1/250 - 1/2000), FC (1/250 - 1/2000)
Expiration Date	12 months from date of receipt.