

Product Datasheet

Goat Anti-Rat IgG (H&L) - Biotin (orb1878924)

| | |
|----------------------------|--|
| Description | Goat Polyclonal Secondary Antibody to Rat IgG (H&L) Biotin |
| Reactivity | Rat |
| Conjugation | Biotin |
| Tested Applications | ELISA, IH, WB |
| Immunogen | Rat IgG |
| Form/Appearance | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Source | Goat |
| Dilution Range | E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |