

Product Datasheet

Goat Anti-Rat IgG (H&L) (orb1878926)

Description	Goat Polyclonal Secondary Antibody to Rat IgG (H&L)
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	Rat IgG
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.