



Product Datasheet

Anti-RSAD2 Antibody (orb1878953)

Description	Rabbit polyclonal antibody to RSAD2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human RSAD2. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8CBB9 , O70600 , Q8WXG1
Entrez	91543 , 65190 , 58185
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of RSAD2 expression...



Immunofluorescent analysis of RSAD2 staining...