


Product Datasheet

Anti-HNRNPA3 Antibody (orb1880029)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to HNRNPA3 |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human HNRNPA3. The exact sequence is proprietary. |
| Form/Appearance | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Source | Rabbit |
| Uniprot ID | P51991 |
| Entrez | 220988 |
| Dilution Range | WB (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HNRNPA3 express...