


Product Datasheet

Anti-CNOT4 Antibody (orb1880392)


Description	Rabbit polyclonal antibody to CNOT4
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	Recombinant fusion protein of human CNOT4. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	095628, Q8BT14
Entrez	53621, 4850
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), FC (1/10 - 1/30)
Expiration Date	12 months from date of receipt.



Western blot analysis of CNOT4 expression...



Immunohistochemical analysis of CNOT4 staining...



Flow cytometric analysis of U2OS cells U...