



## Product Datasheet

### Anti-PER2 Antibody (orb1881019)

<b>Description</b>	Rabbit polyclonal antibody to PER2
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IH, WB
<b>Immunogen</b>	Recombinant fusion protein of human PER2. The exact sequence is proprietary.
<b>Form/Appearance</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">O54943</a> , <a href="#">Q9Z301</a> , <a href="#">O15055</a>
<b>Entrez</b>	<a href="#">63840</a> , <a href="#">18627</a> , <a href="#">8864</a>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IH (1/50 - 1/100)
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of PER2 expression...



Immunohistochemical analysis of PER2 sta...