


Product Datasheet

Gamma-Aminobutyric Acid ELISA Kit (orb1881031)

Description	Competitive Inhibition Enzyme Immunoassay for quantitative detection of Gamma-Aminobutyric Acid conc...
Reactivity	All
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Sample Types	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
Sensitivity	11 pg/ml
Expiration Date	Please enquire.



Representative
Standard
Curve