


Product Datasheet

Dopamine ELISA Kit (orb1881032)

Description	Competitive Inhibition Enzyme Immunoassay for quantitative detection of Dopamine concentrations
Reactivity	All
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Sample Types	Serum, Plasma, Tissue Homogenates or Other Biological Fluids
Sensitivity	7 pg/ml
Expiration Date	Please enquire.



Representative
Standard
Curve