


Product Datasheet

5-Hydroxytryptamine ELISA Kit (orb1881033)

Description	Competitive Inhibition Enzyme Immunoassay for quantitative detection of 5-Hydroxytryptamine concentr...
Reactivity	All
Range	15.6 ng/ml - 1000 ng/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Sample Types	Serum, Plasma or Other Biological Fluids
Sensitivity	4 ng/ml
Expiration Date	Please enquire.



Representative
Standard
Curve