



Product Datasheet

3-Phosphoglycerate Kinase Microplate Assay Kit (orb1881041)

Description	Detection and Quantification of 3-Phosphoglycerate Kinase Activity.
Reactivity	All
Range	3 µmol/L - 300 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-83-6
Expiration Date	Please enquire.



Representative Standard Curve



Positive Control reaction in 96-well pla...