

## Product Datasheet

### Monkey IL-6 (Interleukin 6) ELISA Kit (orb1881450)

<b>Description</b>	Monkey IL-6 (Interleukin 6) ELISA Kit
<b>Reactivity</b>	Monkey
<b>Range</b>	0.313-20pg/ml
<b>Target</b>	IL-6
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, cell culture supernatant and other biological samples
<b>Uniprot ID</b>	<b>P79341</b>
<b>Sensitivity</b>	0.188pg/ml
<b>Expiration Date</b>	Please enquire.