


Product Datasheet

Rat H3cit-NE (neutrophil extracellular trap) ELISA Kit (orb1881472)

Description	Rat H3cit-NE (neutrophil extracellular trap) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	H3cit-NE
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant and other biological samples
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)