

Product Datasheet

Mouse TNFRSF12A (Tumor necrosis factor receptor superfamily member 12A) ELISA Kit (orb1881531)

Description	Mouse TNFRSF12A (Tumor necrosis factor receptor superfamily member 12A) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	TNFRSF12A
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant and other biological samples
Uniprot ID	Q9CR75
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.