

Product Datasheet

Mouse IL-22BP (Interleukin-22-binding protein) ELISA Kit (orb1881536)

Description	Mouse IL-22BP (Interleukin-22-binding protein) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/ml
Target	IL-22BP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant and other biological samples
Uniprot ID	Q80XF5
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.