

Product Datasheet

Mouse IGF-2 (Insulin-like growth factor 2) ELISA Kit (orb1881558)

Description	Mouse IGF-2 (Insulin-like growth factor 2) ELISA Kit
Reactivity	Mouse
Range	1.563-100ng/ml
Target	IGF-2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant and other biological samples
Uniprot ID	P09535
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.