

## Product Datasheet

# Human DBP (Vitamin D Binding Protein) ELISA Kit (orb1881616)

<b>Description</b>	Human DBP (Vitamin D Binding Protein) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.391-25ng/ml
<b>Target</b>	DBP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P02774</b>
<b>Sensitivity</b>	0.234ng/ml
<b>Expiration Date</b>	Please enquire.