
Product Datasheet

Recombinant Mouse Amyloid beta A4 protein (App) (orb1881783)

Description

Recombinant Mouse Amyloid beta A4 protein (App)

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

C-terminal 6xHis-tagged

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Protein Sequence

LEVPTDGNAGLLAEPQIAMFCGKLNMHMNVQNGKWESDPSGTKTCIGTKEGIL
 QYCQEVYPELQITNVVEANQPVTIQNWCKRGRKQCKTHTHIVIPYRCLVGEFVSD
 ALLVPDKCKFLHQERMDVCETHLHWHTVAKETCSEKSTNLHDYGMLLPCGIDKF
 RGVEFVCCPLAEESDSVDSADAEEEDSDVWWGGADTDYADGGEDKVVEVAE
 EEEVADVEEEEADDEDVEDGDEVEEEAEEPYEATERTTSTATTTTTTTSVEE
 VVREVCSEAETGPCRAMISRWFYFDVTEGKCVFFYGGCGGNRNNFDTEEYCM
 AVCGSVSTQSLKTTSEPLPQDPDKLPTTAASTPDAVDKYLETPGDENEHAHFQ
 KAKERLEAKHRERMSQVMREWEAERQAKNLPKADKKAVIQHFQEKVESLEQE
 AANERQQLVETHMARVEAMLNDRRLALENYITALQAVPPRPHHVFNMLKKYVR
 AEQKDRQHTLKHFEHVRMVDPKKAAQIRSQVMTHLRVIYERMNQSLSLLYNVPA
 VAEIQDEVDPELLQKEQNYSDDLANMISEPRISYGNLMPSLTETKTTVELLPV
 NGEFSLDDLQPWHPFGVDSVPANTENEVEPVDARPAADRGLTTRPGSGLTNIKT
 EEISEVKMDAEFGHDSGFVVRHQKLVFFAEDVGSNKGAIIGLMVGGVVIATVIVIT
 LVMLKKKQYTSIHGGVVEVDAAVTPEERHLSKMQQNGYENPTYKFFEQMQN

Purity

Greater than 95% as determined by SDS-PAGE.

Source

E.coli

MW

91.8 kDa

Uniprot ID
[P12023](#)
Expiration Date

6 months from date of receipt.