

Product Datasheet

Recombinant Human Myosin regulatory light chain 12B(MYL12B) (Active) (orb1881866)

Description	Recombinant Human Myosin regulatory light chain 12B(MYL12B) (Active)
Endotoxins	Less than 1.0 EU/ug as determined by LAL method.
Preservatives	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	MSSKKAKTKTKKRPQRATSNVFMFDQSQIQEFKEAFNMIDQNRDGFIDKEDL HDMLASLGKNPTDAYLDAMMNEAPGPINFTMFLTMFGKLNKGTDPEDVIRNAFA CFDEEATGTIQEDYLRELLTTMGDRFTDEEVDELYREAPIDKKGNFNIEFTRILK HGAKDKDD
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	26.7 kDa
Uniprot ID	O14950
Expiration Date	6 months from date of receipt.