

Product Datasheet

Anti-Collagen I/COL1A1 Antibody (orb1882073)

Description	Anti-Collagen I/COL1A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IHC-Fr, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of mouse Collagen I, identical to the related rat sequence, and different from the related human sequence by two amino acids.
Preservatives	5mg rAlbumin, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₆
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water = 500 µg/ml.
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Dilution Range
Expiration Date	12 months from date of receipt.