

Product Datasheet

Anti-Aconitase 2 Antibody Picoband (monoclonal, 4C12D1) (orb1882102)

Description	Anti-Aconitase 2 Antibody Picoband (monoclonal, 4C12D1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Aconitase 2, identical to the related mouse and rat sequences.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.