

Product Datasheet

Anti-MDH2 Antibody Picoband (monoclonal, 5D8C1) (orb1882105)

Description	Anti-MDH2 Antibody Picoband (monoclonal, 5D8C1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	E.coli-derived human MDH2 recombinant protein (Position: A9-L223). Human MDH2 shares 97.7% and 98.1% amino acid (aa) sequence identity with mouse and rat MDH2, respectively.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.