

Product Datasheet

Anti-Claudin 3/CLDN3 Antibody Picoband (monoclonal, 4C7D2) (orb1882108)

Description	Anti-Claudin 3/CLDN3 Antibody Picoband (monoclonal, 4C7D2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	E.coli-derived human Claudin 3/CLDN3 recombinant protein (Position: A101-V220).
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.