

Product Datasheet

Anti-Ataxin 1 Antibody Picoband (monoclonal, 2B13G8) (orb1882116)

Description	Anti-Ataxin 1 Antibody Picoband (monoclonal, 2B13G8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Ataxin 1, different from the related mouse and rat sequences by one amino acid.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.