

Product Datasheet

Anti-Drebrin/DBN1 Antibody Picoband (monoclonal, 4F6E7) (orb1882128)

| | |
|----------------------------|---|
| Description | Anti-Drebrin/DBN1 Antibody Picoband (monoclonal, 4F6E7) |
| Species/Host | Mouse |
| Reactivity | Human, Monkey, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | FC, IHC, WB |
| Immunogen | E.coli-derived human Drebrin/DBN1 recombinant protein (Position: H9-D649). |
| Preservatives | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ . |
| Form/Appearance | Lyophilized |
| Concentration | Adding 0.2 ml of distilled water =f 500 µg/ml. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Dilution Range | „ |
| Expiration Date | 12 months from date of receipt. |