

Product Datasheet

Anti-Sorbitol Dehydrogenase/SORD Antibody Picoband (monoclonal, 12B10G2) (orb1882153)

Description	Anti-Sorbitol Dehydrogenase/SORD Antibody Picoband (monoclonal, 12B10G2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	E.coli-derived human Sorbitol Dehydrogenase/SORD recombinant protein (Position: N8-P357).
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.