

Product Datasheet

Anti-HGS Antibody Picoband (monoclonal, 6C2E1) (orb1882156)

| | |
|----------------------------|--|
| Description | Anti-HGS Antibody Picoband (monoclonal, 6C2E1) |
| Species/Host | Mouse |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | FC, IHC, WB |
| Immunogen | E.coli-derived human HGS recombinant protein (Position: R3-D777). |
| Preservatives | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ . |
| Form/Appearance | Lyophilized |
| Concentration | Adding 0.2 ml of distilled water =f 500 µg/ml. |
| Storage | -20°C |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Dilution Range | „ |
| Expiration Date | 12 months from date of receipt. |