

## Product Datasheet

# Anti-RALBP1 Antibody Picoband (monoclonal, 2F9D3) (orb1882165)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Anti-RALBP1 Antibody Picoband (monoclonal, 2F9D3)  |
| <b>Species/Host</b>        | Mouse  |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | FC, ICC, IF, IHC, WB   |
| <b>Immunogen</b>           | E. coli-derived human RALBP1 recombinant protein (Position: K239-Q506).  |
| <b>Preservatives</b>       | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> . |
| <b>Form/Appearance</b>     | Lyophilized  |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water =f 500 µg/ml.   |
| <b>Storage</b>             | -20°C  |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG2b  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Dilution Range</b>      | ....   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |