

Product Datasheet

Anti-p95 NBS1 Antibody Picoband (monoclonal, 5D7A1) (orb1882189)

Description	Anti-p95 NBS1 Antibody Picoband (monoclonal, 5D7A1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human p95 NBS1, different from the related mouse sequence by three amino acids, and from the related rat sequence by five amino acids.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.