

Product Datasheet

Anti-PDE6 beta/PDE6B Antibody Picoband (monoclonal, 8I2D7) (orb1882193)

Description	Anti-PDE6 beta/PDE6B Antibody Picoband (monoclonal, 8I2D7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human PDE6 beta/PDE6B recombinant protein (Position: K25-Q237).
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range
Expiration Date	12 months from date of receipt.