


Product Datasheet

ALOX5AP antibody (orb18822)

Description	Goat polyclonal antibody to ALOX5AP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	FLAP / ALOX5AP (C Term)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: This product has been successfully used in transfected HeLa: (Maier et al, Biochem Pharmacol. 2008 Oct 1;76(7):862-72.; PubMed: 18692027). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KTISTTISPLLLIP
Clonality	Polyclonal
MW	18.2; 24.2
Entrez	241
Dilution Range	ELISA: 1:16000, IHC-P: 6 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Sp...



Immunohistochemical staining of Human Sp...