

## Product Datasheet

### Anti-Ki67 Antibody Picoband (monoclonal, 5E12) (orb1882218)

<b>Description</b>	Anti-Ki67 Antibody Picoband (monoclonal, 5E12)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IF, IHC
<b>Immunogen</b>	E. coli-derived human Ki67 recombinant protein (Position: K2860-I3256).
<b>Preservatives</b>	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water =f 500 µg/ml.
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	...
<b>Expiration Date</b>	12 months from date of receipt.