

## Product Datasheet

# Anti-Ku70 Antibody Picoband (monoclonal, 9B6) (orb1882228)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Anti-Ku70 Antibody Picoband<br>(monoclonal, 9B6)  |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | FC, ICC, IF, IHC, WB  |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence at C-terminus of human Ku70, different from the related mouse sequence by one amino acid. |
| <b>Preservatives</b>       | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> .                          |
| <b>Form/Appearance</b>     | Lyophilized   |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water =f 500 µg/ml.  |
| <b>Storage</b>             | -20°C   |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG2b   |
| <b>Clonality</b>           | Monoclonal  |
| <b>Dilution Range</b>      | ....  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |