

## Product Datasheet

# Anti-Fibrinogen beta chain/FGB Antibody Picoband (monoclonal, 6D12) (orb1882256)

<b>Description</b>	Anti-Fibrinogen beta chain/FGB Antibody Picoband (monoclonal, 6D12)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	A synthetic peptide corresponding to a sequence in the middle region of human FGB, different from the related mouse sequence by three amino acids, and from the related rat sequence by five amino acids.
<b>Preservatives</b>	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> .
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water =f 500 µg/ml.
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	,
<b>Expiration Date</b>	12 months from date of receipt.