

## Product Datasheet

# Anti-Liver Carboxylesterase 1/CES1 Antibody Picoband (monoclonal, 3F10) (orb1882265)

<b>Description</b>	Anti-Liver Carboxylesterase 1/CES1 Antibody Picoband (monoclonal, 3F10)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IF, IHC, WB
<b>Immunogen</b>	E. coli-derived human Liver Carboxylesterase 1/CES1 recombinant protein (Position: E99-A206).
<b>Preservatives</b>	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Na <sub>3</sub> .
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water =f 500 µg/ml.
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	....
<b>Expiration Date</b>	12 months from date of receipt.