

Product Datasheet

Anti-TIF1 gamma/TRIM33 Antibody Picoband (monoclonal, 8I8) (orb1882276)

Description	Anti-TIF1 gamma/TRIM33 Antibody Picoband (monoclonal, 8I8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human TIF1 gamma recombinant protein (Position: M1001-K1127). Human TIF1 gamma shares 96.1% amino acid (aa) sequence identity with mouse TIF1 gamma.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.