

## Product Datasheet

# Anti-PP2A-alpha/PPP2CA Antibody Picoband (monoclonal, 3B6) (orb1882277)

<b>Description</b>	Anti-PP2A-alpha/PPP2CA Antibody Picoband (monoclonal, 3B6)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IF, IHC, WB
<b>Immunogen</b>	E.coli-derived human PP2A-alpha recombinant protein (Position: M1-L309). Human PP2A-alpha shares 100% amino acid (aa) sequence identity with both mouse and rat PP2A-alpha.
<b>Preservatives</b>	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Na <sub>3</sub> .
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water =f 500 µg/ml.
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	,,,,
<b>Expiration Date</b>	12 months from date of receipt.