

Product Datasheet

Anti-Bid Antibody Picoband (monoclonal, 10F4) (orb1882278)

Description	Anti-Bid Antibody Picoband (monoclonal, 10F4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human Bid recombinant protein (Position: M1-D195). Human Bid shares 64% and 61% amino acid (aa) sequences identity with mouse and rat Bid, respectively.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.