

Product Datasheet

Anti-CD146/MCAM Antibody Picoband (monoclonal, 4C12) (orb1882285)

Description	Anti-CD146/MCAM Antibody Picoband (monoclonal, 4C12)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-Fr, WB
Immunogen	E.coli-derived human CD146 recombinant protein (Position: H59-A401). Human CD146 shares 73% amino acid (aa) sequence identity with both mouse and rat CD146.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.