


## Product Datasheet

### Pescadillo antibody (orb18841)

<b>Description</b>	Goat polyclonal to Pescadillo
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Bovine, Human, Mouse, Porcine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Target</b>	Pescadillo
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 70kDa band observed in HeLa lysate (calculated MW of 68kDa according to NP_055118). Recommended concentration: 1-3 µg/ml.
<b>Protein Sequence</b>	GGLEKKKYERGSA
<b>Clonality</b>	Polyclonal
<b>MW</b>	68
<b>Entrez</b>	<a href="#">23481</a>
<b>Dilution Range</b>	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 2.5-3.8µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa lysate usi...



Immunohistochemical staining of Human Te...