

Product Datasheet

PACSIN1 antibody (orb18842)

Description

Goat polyclonal antibody to PACSIN1

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Porcine

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

PACSIN1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Human Brain (hippocampus) lysates (calculated MW of 51.1kDa according to NP_065855.1). Recommended concentration: 0.03-0.1 µg/ml. An additional band of unknown identity was also consistently observed at 37kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

SSSYDEASLAPEET

Clonality

Polyclonal

MW

51.1

Entrez


[29993](#)

Dilution Range


ELISA: 1:32000, WB: 0.03-0.1 µg/ml, IHC-P: 2.5 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of Hippocampus lys...



Immunohistochemical staining of Human Ce...



Immunohistochemical staining of Human Co...