
Product Datasheet

Histone H3 (phospho-Thr45) antibody (FITC) (orb188482)

Description

Rabbit polyclonal antibody against

conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Thr45

Target

H3-4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

17kDa

Uniprot ID[Q16695](#)**Dilution Range**

IF=1:100-500,ICC=1:100-500

Expiration Date

12 months from date of receipt.