

Product Datasheet

Influenza A H1N1 Neuraminidase antibody (FITC) (orb188645)

Description

Rabbit polyclonal antibody against

conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide
derived from Influenza A H1N1
Neuraminidase

Target

H1N1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid
repeated freeze/thaw cycles. The
lyophilized antibody is stable at
room temperature for at least one
month and for greater than a year
when kept at -20°C. When
reconstituted in sterile pH 7.4 0.01M
PBS or diluent of antibody the
antibody is stable for at least two
weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

52kDa

Uniprot ID
A8C8J7
Dilution Range

IF=1:100-500, ICC=1:100-500

Expiration Date

12 months from date of receipt.