


Product Datasheet

STAP2 antibody (orb18881)

Description	Goat polyclonal antibody to STAP2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	STAP2 / BKS
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in human heart lysates (calculated MW of 49.6kDa according to NP_060190.2). Recommended concentration: 0.05-0.3 µg/ml.
Protein Sequence	ELQKKLEKRRALEH
Clonality	Polyclonal
MW	49.6; 44.9
Entrez	55620
Dilution Range	ELISA: 1:128000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart lys...



Immunohistochemical staining of Human Co...