


Product Datasheet

NKH477 (orb1888129)

Description	NKH477 (Colforsin dapropate hydrochloride) is a novel water-soluble forskolin derivative that improv...
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles
Note	For research use only
Formula	C ₂₇ H ₄₄ ClNO ₈
MW	546.09
CAS Number	138605-00-2
SMILES	<chem>O[C@@]12[C@@]3([C@@](C(C)(CC[C@H]3O)C)([H])[C@H](OC(CCN(C)C)=O)[C@H](OC(C)=O)[C@]1(O[C@](C=C)(C)CC2=O)C)C.Cl</chem>
Expiration Date	12 months from date of receipt.



NKH477
Chemical
Structure