

Product Datasheet

PRDM4 antibody (orb18914)

Description	Goat polyclonal antibody to PRDM4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRDM4 / PFM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 70kDa band observed in human and mouse brain lysates (calculated MW of 87.9kDa according to NP_036538). Recommended concentration: 0.3-2 µg/ml.
Protein Sequence	CSAVYSADESLSAHK
Clonality	Polyclonal
MW	87.9
Entrez	11108
Dilution Range	ELISA: 1:64000, WB: 0.3-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18914
(0.3ug/ml)
staining of
human br...