
Product Datasheet

SYNGR3 antibody (orb18925)

Description

Goat polyclonal antibody to SYNGR3

Species/Host

Goat

Reactivity

Canine, Human, Rat

Conjugation

Unconjugated

Tested Applications

ELISA

Target

Synaptogyrin 3

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave an approx 50kDa band in Human and Mouse Brain lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 24.6kDa according to NP_004200. The 50kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

QRTAPGPATTQAGD

Clonality

Polyclonal

MW

24.6

Entrez
9143
Dilution Range

ELISA: 1:32000, IHC-P: 10ug/ml

Expiration Date

12 months from date of receipt.