


Product Datasheet

SYNGR4 antibody (orb18942)

Description	Goat polyclonal antibody to SYNGR4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Synaptogyrin 4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 23kDa band observed in Human Brain lysates (calculated MW of 25.8kDa according to NP_036583). Recommended concentration: 0.01-0.1 µg/ml. A minor band of unknown identity was also consistently observed at 150kDa. This band was successfully blocked by incubation with the immunising peptide.
Protein Sequence	HIPKSLQELANS
Clonality	Polyclonal
MW	25.8
Entrez	23546
Dilution Range	ELISA: 1:128000, WB: 0.01-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb18942
(0.01µg/ml)
staining of
Human B...