

Product Datasheet

NFkB p65 (phospho-Ser281) antibody (FITC) (orb189970)

Description	NFkB p65 (phospho-Ser281) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NFkB p65 around the phosphorylation site of Ser281
Target	phospho-NFkB p65 (Ser281)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04206
Expiration Date	12 months from date of receipt.